

SupraScooter® MicroLite

THE WORLDS LIGHTEST MOBILITY SCOOTER

A stunning design continues the tradition of innovation from our design team and was inspired by the need to provide safe and reliable portable transport far beyond what was thought possible. Achieving a staggeringly low weight means little if the performance is not equally impressive. Weighing only 18kgs the MicroLite ticks all the boxes for easy transportation, amazing indoor manoeuvrability and incredible outdoor performance.

TEST DRIVE A MICROLITE AND YOUR PERCEPTION OF WHAT A TRULY PORTABLE MOBILITY SCOOTER IS CAPABLE OF WILL CHANGE FOREVER



Australian Concept Design

Specifications*

Length	930mm
Width	500mm
Max Load (Level Ground)	100kg
Max Speed	6kph
Est. Travel Range	8 - 10km
Tyre Sizes (Solid)	175 x 50mm
Battery (Lithium)	24V 11.6Ah
Motor	230w
Total Weight	18kg
Weight Disassembled	11.8kg
Front Basket	Optional

*Specifications subject to change without notice

trekmobility.com.au