

# TREK BRONCO

## Experience the Power of the Bronco

Unleash your adventurous spirit with the Bronco, an all terrain scooter engineered for performance and durability. Whether you're navigating outdoor landscapes or cruising through urban environments, the Bronco combines cutting edge technology with stunning design to deliver a high class driving experience. Get ready to conquer any terrain with confidence and style.



### Specifications\*

Motor	YC 1400W(Taiwan)
Controller	PG Drive200A
Battery	75A to 100A
Charger	8A
Max Speed	10km/h
Weight Capacity	250kg
Climbing	12°
Turning Radius	1815mm
Size	1600*805*1435mm
Tire Size	13*6.50-6 Tubeless
Weight(WO/batteries)	120kg

\*subject to change without notice