

# TREK BRAVO

## POWER, AGILITY AND A SUPERB DRIVING EXPERIENCE

ENGINEERED FOR HIGH PERFORMANCE AND RESPONSIVE HANDLING

Our best selling scooter for 5 years in a row, the Bravo offers an unrivalled blend of comfort, refinement and class leading features at a highly competitive price. The powerful 4 pole motor and transaxle have been carefully matched to provide stunning hill climbing performance and has equally impressive low speed control and maneuverability when negotiating busy shopping strips and markets.



### Specifications\*

Length	1400mm
Width	686mm
Max Load (On Flat Ground)	185kg
Max Speed programmable up to	15km/h
Travel range up to...	45km
Suspension front & rear	YES
Tyre Sizes	330mm
Battery	55ah
Charger	6a
Motor (Taiwanese 4 pole)	
Controller PG S-Drive	140a

### Features

- \* Seat Belt
- \* Free Rear Basket
- \* High/Low Speed Switch
- \* Reclining Back Rest
- \* LED Lighting & Indicators
- \* Front and Rear Suspension
- \* Seat Suspension

\*subject to change without notice

[trekmobility.com.au](http://trekmobility.com.au)