

TREK Evolution 4

COMBINING ELEGANT SIMPLICITY WITH AGILITY
AND UNMATCHED RELIABILITY

UNRIVALLED COMFORT REFINEMENT AND PERFORMANCE

Possibly the best overall luxury scooter value in its class. The Evolution 4 is ideally suited to users who require a feature rich luxury product at an attractive price.

Designed by Korea's leading scooter manufacturer, features include a powerful 4 pole motor, front and rear suspension, LED lighting and fully adjustable high back seat.



Specifications*

Length	1345MM
Width	695MM
Max Load	160KG
Max Speed programmable up to	12KM/H
Travel range up to...	40KM
Suspension front & rear	YES
Tyre Sizes	300MM
Battery	45AH
Charger	6A
Motor (Taiwanese 4 pole)	
Controller PG S-Drive	120A

*subject to change without notice

Features

- * Seat Belt
- * Free Rear Basket
- * High/Low Speed Switch
- * Reclining Back Rest
- * LED Lighting & Indicators
- * Front and Rear Suspension
- * Seat Suspension