

# SupaScooter® SupaLite

## TURNING VISION INTO REALITY

WHEN INNOVATION AND INSPIRATION COLLIDED THE SUPALITE WAS BORN

The SupaLite is more than a product, it's about innovation and performance that inspires and fulfils a need. Look around, there is nothing like it. It's not the result of marketing men, committees or audience clinics. A well designed product is one you fall in love with, it has to provoke a real reaction. Our design team have a history of innovation and creating products that customers want to buy. The SupaLite is a product you want to engage with now and in the future.

THE SUPALITE WAS BORN NOT OF THE OLD SCHOOL OR THE NEW SCHOOL, BUT OF IT'S OWN SCHOOL. TEST RIDE AND EXPLORE THE POSSIBILITIES, YOU WILL BE AMAZED.



### Specifications\*

Length	970mm
Width	570mm
Max Load (Level Ground )	120kg
Max Speed	6.5kp/h
Est. Travel Range	10 - 13km
Pneumatic Tyres	200 x 65mm
Battery (Lithium)	25.2V 14Ah
Motor	260w
Total Weight	20.4kg
Weight (Disassembled)	13.6kg
Front Basket	Optional
Arm Rests	Optional

\*Specifications subject to change without notice

[trekmobility.com.au](http://trekmobility.com.au)