

TREK BRAVO

POWER, AGILITY AND A SUPERB DRIVING EXPERIENCE

ENGINEERED FOR HIGH PERFORMANCE AND RESPONSIVE HANDLING

Our best selling scooter for 5 years in a row, the Bravo offers an unrivalled blend of comfort, refinement and class leading features at a highly competitive price. The powerful 4 pole motor and transaxle have been carefully matched to provide stunning hill climbing performance and has equally impressive low speed control and maneuverability when negotiating busy shopping strips and markets.



Specifications*

Length	1400mm
Width	686mm
Max Load (On Flat Ground)	185kg
Max Speed programmable up to	15km/h
Travel range up to...	45km
Suspension front & rear	YES
Tyre Sizes	330mm
Battery	55ah
Charger	6a
Motor (Taiwanese 4 pole)	
Controller PG S-Drive	140a

Features

- * Seat Belt
- * Free Rear Basket
- * High/Low Speed Switch
- * Reclining Back Rest
- * LED Lighting & Indicators
- * Front and Rear Suspension
- * Seat Suspension

*subject to change without notice

trekmobility.com.au