

TREK Safari

Mobility in Total Comfort

The Trek Safari blends rugged capability with a lightweight, portable design that disassembles in seconds for easy transport. Built for all-day support, it features a high-back seat with an adjustable backrest, integrated seat suspension, and a 360-degree swivel for effortless access at home or on the go.



Specifications*

- Length 910mm
- Width 610mm
- Height 865mm
- Swivel Seat
- Weight Capacity 125kg
- Max Speed 9km
- Turning Radius 730mm
- Motor Unit Weight 20.2kg
- Seat Weight 10.6kg
- Battery 24V 15AH

Co-designed by
Trek Mobility & Eurogreen International
Manufactured in Taiwan



*subject to change without notice

trekmobility.com.au