

TREK rhino

Embrace the journey with our newest addition

Designed as a budget-friendly alternative to the Trek Evolution 4, the Rhino packs all the essential features you need for comfortable and reliable mobility. With its sturdy build, smooth maneuverability, and long-lasting battery life, the Rhino ensures you can navigate through your day with ease and confidence.



Specifications*

Colours	Red & Grey
Length	1345MM
Width	695MM
Max Load	150KG
Max Speed	10KM/H
Travel range up to...	40KM
Tyre Sizes	300MM
Battery	45AH
Charger	6A
Motor (Taiwanese 4 pole)	500W
Controller PG S-Drive	90A

*subject to change without notice